

Part No.


111 8050H500S

Description

Disposable LiquiThief – 500mm

Material of Construction:

Tip	HDPE (virgin, FDA approved – CFR 177.1520)
Outer Tube	HDPE (virgin, FDA approved – CFR 177.1520)
Inner Rod	HDPE (virgin, FDA approved – CFR 177.1520)
Clutch	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Housing	HDPE (virgin, FDA approved – CFR 177.1520)
Piston Seal	Silicone Rubber (virgin, FDA approved – CFR 177.2600)
Handle	HDPE (virgin, FDA approved – CFR 177.1520)



Method of Construction:

Ultrasonic welding

Nominal Length

Fully closed:	523mm
Fully extended:	881mm

Diameter of Outer Tube:

21mm

Diameter of hole in tip:

1.5mm

Nominal Weight of Sampler:

83g

Assembly & Packing Environment:

Class 100,000 Medical Cleanroom

Individually Bagged?

Yes (heat sealed PE bag)

Method of Sterilisation:

Gamma Irradiation (25-45kGy)

Number of Samplers per Box:

20

Recommended Storage Conditions:

Dry and ambient temperature